ZDOCK195 Rev. 20 - FBS Kits with Sampling Wand

Language		Page
EN		2-7
DE		8-14
DA		15-20
ES		21-26
FI		27-32
FR		33-39
NL		40-46
SV		47-52

CE
1639

ROCKET MEDICAL PLC
Sedling Road, Washington,
NE38 9BZ, UK

EC REP
ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL LLC
50 Corporate Park Drive #890.
East Pembroke. MA. 02359.
USA

Australian Sponsor:
ROCKET MEDICAL PTY LTD
Suite 209, 1 Katherine Street,
Chatswood, NSW 2067 Australia

The Rocket device and the Rocket Medical logo are registered trademarks of Rocket Medical Plc.

Product Name: Rocket Medical Fetal Blood Sampling (FBS) Kits


Procedure Pack Product Codes:

R57018-00-SW - FBS Kit - 5 x 55µL Capillary Tubes

R57024-00-SW - FBS Kit - 5 x 85µL Capillary Tubes

R57026-00-SW - FBS Kit - 5 x 100µL Capillary Tubes

Device Image:

Device Contents:

DEVICE REF		DEVICE NAME	QTY	LEGAL MANUFACTURER DETAILS
R57018-00-SW	ZPACK011	FBS Sampling Wand for 55µl Capillary Tubes	1	**ROCKET MEDICAL PLC** Sedling Road, Washington, England, NE38 9BZ www.rocketmedical.com
		Standard FBS Amnioscope	1	
		Amnilume™	1	
		Capillary Tube Holder	1	
		Paraffin Tube	1	
		Capillary Tube Caps	10	
		Gauze Swabs	5	
	ZTUBE095	55µl Capillary Tube	5	**RADIOMETER MEDICAL APS**, Åkandevej 21, 2700 København, Denmark radiometer@radiometer.dk
	ZSPNG007	Sponge Holding Forceps	1	**WILLIAMS MEDICAL SUPPLIES LTD**, Craiglas House, Maerdy Industrial Estate North, Rhymney, NP22 5PY, UK customercare@wms.co.uk
	ZSWAB008	Swab Sticks	6	**APPLIMED SA**, Rte de Pra de Plan 1, 1618 Châtel-Saint-Denis, Switzerland info@applimed.ch
R57024-00-SW	ZPACK012	FBS Sampling Wand for 85µl Capillary Tubes	1	**ROCKET MEDICAL PLC** Sedling Road, Washington, England, NE38 9BZ www.rocketmedical.com
		Standard FBS Amnioscope	1	
		Amnilume™	1	
		Capillary Tube Holder	1	
		Paraffin Tube	1	
		Capillary Tube Caps	10	
		Gauze Swabs	5	

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

 Copyright © 2023 All rights reserved.

	ZTUBE094	85µl Capillary Tube	5	**RADIOMETER MEDICAL APS**, Åkandevej 21, 2700 København, Denmark radiometer@radiometer.dk
	ZSPNG007	Sponge Holding Forceps	1	**WILLIAMS MEDICAL SUPPLIES LTD**, Craiglas House, Maerdy Industrial Estate North, Rhymney, NP22 5PY, UK customercare@wms.co.uk
	ZSWAB008	Swab Sticks	6	**APPLIMED SA**, Rte de Pra de Plan 1, 1618 Châtel-Saint-Denis, Switzerland info@applimed.ch
R57026-00-SW	ZPACK013	FBS Sampling Wand for 100µl Capillary Tubes	1	**ROCKET MEDICAL PLC** Sedling Road, Washington, England, NE38 9BZ www.rocketmedical.com
		Standard FBS Amnioscope	1	
		Amnilume™	1	
		Capillary Tube Holder	1	
		Paraffin Tube	1	
		Capillary Tube Caps	10	
		Gauze Swabs	5	
	ZTUBE113	100µl Capillary Tube	5	**VITREX MEDICAL A/S**, Vasekær 6, 2730 Herlev, Denmark vitrex@vitrex.dk
	ZSPNG007	Sponge Holding Forceps	1	**WILLIAMS MEDICAL SUPPLIES LTD**, Craiglas House, Maerdy Industrial Estate North, Rhymney, NP22 5PY, UK customercare@wms.co.uk
	ZSWAB008	Swab Sticks	6	**APPLIMED SA**, Rte de Pra de Plan 1, 1618 Châtel-Saint-Denis, Switzerland info@applimed.ch

Intended Use:

The sampling wand is intended to be used to make a small incision in the fetal scalp and to hold the capillary tube in order to collect fetal blood sample(s) for subsequent analysis.

The Rocket Medical Amnioscope and Amnilume™ are intended to provide access to and illumination of the fetal scalp during fetal blood sampling.

The Rocket Medical capillary tube holder is intended to be used to hold the capillary tube(s) in order to collect fetal blood sample(s) for subsequent analysis.

The Rocket Medical paraffin is intended to be used by application to the fetal scalp to aid blood pooling/droplet formation for collection for subsequent analysis.

The Rocket Medical capillary tube cap is intended to be used to protect the fetal blood sample during transfer to the third-party analyser.

The Rocket Medical gauze swabs are intended to be used to clean the fetal scalp.

Indications:

The Rocket Medical devices in this procedure pack are indicated for use for patients requiring fetal blood analysis (for pH/blood gas levels) during labour. Dilation of at least 3cm is required, with ruptured membranes.

Contraindications:

Fetal blood sampling is contraindicated:

- when there is acute fetal compromise (for example, fetal bradycardia of >3 minutes)
- when there is an acute event such as cord prolapse, uterine rupture or suspected placental abruption
- where the full clinical picture suggests that the birth should be expedited, for example, in cases of maternal sepsis or unstable pre-eclampsia

as this would delay an urgent requirement to deliver the fetus, leading to deterioration in the fetal condition.

FBS is also contraindicated:

- in cases of placenta previa as this may cause serious bleeding.
- in maternal thrombocytopenia as this may lead to increased bleeding, bruising, or abnormal blood clotting.
- when the fetus is premature (less than 34+0 weeks) as this may lead to or worsen existing fetal anaemia.
- where the mother is a confirmed carrier of haemophilia and the fetus is either affected or of unknown status as this may lead to increased bleeding.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

 Copyright © 2023 All rights reserved.

- where there is a risk of maternal-fetal infection transfer (for example, active herpes simplex virus, HIV and hepatitis C). Use caution and refer to your local clinical policies and procedures in cases when the woman is positive for Group B streptococcus, gonorrhoea, hepatitis A, hepatitis B or syphilis.
- in cases of facial presentation as there may be permanent facial scarring or injury to the eye.
- when identification of the presenting part is uncertain as an incision may be made into the fetal genitals, face or fontenelles, leading to permanent injury to the eye, loss of CSF or scaring.

WARNING:	Ensure close contact between amnioscope and fetal scalp is maintained during the sampling procedure to prevent leakage of amniotic fluid into the blood sample, as this could lead to a falsely reassuring result.
WARNING:	Do not use this kit without prior verification that it is suitable for use with the analyser.
WARNING:	There is uncertainty about the usefulness of FBS in women with fever, sepsis or suspected sepsis, since FBS may be falsely reassuring. In these cases, the decision to perform FBS or results obtained from FBS should be discussed with a consultant obstetrician.
WARNING:	This device contains 2 lithium-based primary cells. Do NOT attempt to disassemble the device. Do NOT attempt to recharge or short circuit the batteries. Lithium cells are suitable for incineration or landfill disposal.
CAUTION:	Avoid fetal face, fontanelles and genitalia. Rocket Medical FBS blades have a maximum depth of penetration of 2.1mm
CAUTION:	Do not use Ethyl Chloride Spray – this may adversely affect the internal surface of the plastic amnioscope
CAUTION:	Ensure that both the obturator and Amnilume™ are in position in the amnioscope before inserting into the patient. If the Amnilume™ is not in place, the obturator can protrude beyond the end of the amnioscope, which may cause patient tissue to become trapped between the amnioscope and the obturator.
CAUTION:	Carefully follow the analyser manufacturer's instructions to determine the capillary tube required. Failure to use correct blood volume, select correct analysis mode or the use of clotted or contaminated samples may lead to errors in results or rejection of the whole sample.
CAUTION:	Do NOT use the sample if the blood clots.
CAUTION:	Avoid distributing air through the sample; samples containing air bubbles may be rejected by the analyser.

Benefit Risk: Rocket Medical plc has taken all necessary steps to ensure that residual risks associated with the use of FBS Kits with Sampling Wand are reduced as far as possible through application of existing state of the art techniques in the design and manufacture of these medical devices to ensure safe usage. Rocket Medical plc concludes that the overall medical benefits of FBS Kits with Sampling Wand outweigh the possible risks when used according to the intended use.

Performance Claims and Clinical Benefit: The intended performance claims for the device are:

- To make a small incision in the fetal scalp
- Device is safe to use
- To collect blood suitable for subsequent analysis

The intended clinical benefit is:

- The determination of fetal blood pH level to aid in decision making during labour

Undesirable Side Effects: Fetal side effects include: bleeding, haematoma, infection at incision point.

Description of the Device: The Rocket Medical amnioscope and amnilume connect together to provide access to and illumination of the fetal scalp. The amnioscope - with the amnilume attached - is inserted into the vagina with the plastic obturator in place to provide a smooth leading edge. The obturator (which is considered a device together with the amnioscope) is then removed.

The amnioscope provides access to the fetal head and also assists in preventing ingress of amniotic fluid into the blood sample. Amniotic fluid could change the reading and give a falsely reassuring result.

Once the amnioscope has been inserted, the fetal scalp is prepared for blood collection. The scalp is cleaned using a gauze swab and forceps. Next, a thin layer of paraffin is applied to the fetal scalp (using swab sticks) to aid pooling of the blood sample for collection. The procedure packs contain five tubes from one of three different sizes of capillary tube.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

 Copyright © 2023 All rights reserved.

In preparation for sample collection, a capillary tube is placed into the holder of a Rocket Medical Sampling Wand or a Rocket Medical Capillary Tube Holder, according to the preference of the clinician. The Sampling Wand has a holder for a single tube, whilst the Capillary Tube Holder can hold up to two tubes, one at either end.

The sampling wand is a device that combines the means to puncture the fetal scalp and to hold the capillary tube for blood sample collection by capillary action.

Following scalp preparation, the Sampling Wand is then inserted through the amnioscope. When the fetal scalp has been punctured and a droplet/pool of fetal blood has formed, the capillary tube is presented to the droplet of blood, which enters the tube by capillary action. Capillary tube caps are provided in the procedure pack, which may be used to protect the sample during transfer to the analyser. The capillary tube is inserted into a blood gas analyser, which will provide the pH/blood gas measurement, The measurement is interpreted using clinical guidelines.

One procedure pack may take up to five samples through presentation of multiple capillary tubes to the blood droplet. Once the blood sample(s) has/have been collected, the Rocket Medical FBS devices are removed from the vagina and should be discarded. The blade is not intended to be reused since it may have become blunt and reuse could also present a risk of infection. Reintroduction of the amnioscope may also increase the risk of infection.

User: The intended user is a healthcare professional proficient in the undertaking of fetal blood sampling, working in accordance with local and national guidelines. The device may also be used by a trainee under the supervision of a healthcare professional proficient in the undertaking of FBS.

Patient Population: The intended demographic for this device is pregnant patients during labour.

Environment: To be used in a clinical setting.

For Single Use Only: the device may not have structural integrity to support reuse, resulting in harm to the patient and/or user. Re-use may result in harm to patients and users from cross-contamination and infection.

Sterile: the devices are supplied sterile by Ethylene Oxide (EO) / Gamma Irradiation. Do not re-sterilise; it may not be possible to achieve the required sterility assurance level upon re-sterilisation of a used device. Re-sterilisation may also compromise the structural integrity of the device, leading to device failure.

MRI Safety: This device is not intended to be routinely found in an MRI environment and compatibility has not been verified.

General Information: Rocket Medical supplies a range of FBS Kits suitable for a wide range of blood gas analysers. Please refer to the table below for the specification of the capillary tubes supplied within each kit and check your analyser manual and sampling volume to select the most appropriate kit for your analyser. If more than one tube is suitable, select the smallest volume tube to minimise the blood volume collected.

FBS Kit	Volume µl	Tube Length	OD	ID	Heparin Type	Concentration IU/ml
R57018-00-SW*	55	60 mm +/- 1 mm	2.10 mm + 0.05 mm / - 0.1 mm	1.25 mm + 0.1 mm / -0.05 mm	Na	240
R57024-00-SW*	85	87 mm +/- 1 mm	2.10 mm + 0.05 mm / - 0.10 mm	1.25 mm + 0.10 mm / - 0.05 mm	Na	240
R57026-00-SW	100	85.0mm	2.3 ± 0.05mm	1.20 ± 0.05mm	Ca/Na/Li	200

*Not suitable for lactate measurement on Radiometer analysers as this measurement requires a balanced heparin tube.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

 Copyright © 2023 All rights reserved.

Symbols and safety signs used on the device and labelling:


Symbol or Safety Sign	Meaning	Symbol or Safety Sign	Meaning
	Manufacturer		Keep away from sunlight
EC REP	Authorized representative in the European Community/European Union		Keep dry
Rx ONLY	Caution: Federal law restricts this device to sale by or on the order of a physician.		Do not re-use
	Date of manufacture	MD	Medical device
	Use-by date		Consult the instructions for use
LOT	Batch code		Caution
REF	Catalogue number	STERILIZE	Do not resterilise
STERILE EO	Sterilized using ethylene oxide	LATEX	Does not contain natural rubber latex
Li	Indicates a medical device that contains a Lithium coin cell.		Do not use if package is damaged and consult instructions for use
	Single sterile barrier system with protective packaging outside		Single sterile barrier system
	Distributor		Importer

Operating Instructions:

CHECK: Do not use if the packaging is open or damaged or the devices are damaged. Ensure you have the full contents of the kit prior to use (see 'Device Contents').

USE:

1. Using aseptic technique, remove the kit from outer bag and unfold outer drape to form a sterile field.
2. Select a capillary tube from the tube in the kit and load into the holder on the sampling wand. Locate the end of the tube into the holder at (**1**) and then gently advance into the rear tube grip at (**2**). Ensure the tube is secure before sampling. The kit contains a secondary double-ended capillary tube holder if preferred.


3. Follow local hospital procedure to prepare the perineum and vagina, and complete a vaginal examination to establish the cervical dilatation and station of the fetal head.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

 Copyright © 2023 All rights reserved.

4. Remove the obturator from the amnioscope. Locate the Amnilume™ light module into the slot in the amnioscope and push gently to engage. The light will activate when the light source is correctly fitted. Between samplings, remove the Amnilume™ light source from the endoscope by pushing off with the thumb.
5. Insert amnioscope complete with obturator into the vagina against the fetal presentation ensuring that the cervix is not trapped.
6. Remove obturator. Position the bevelled end by rotating amnioscope to obtain the best possible seal against the fetal presentation. Rotate to finally position the light source module in the upper quadrant.
7. Maintain close contact between amnioscope and scalp to prevent leakage of liquor.
8. Clean puncture area with the green x-ray detectable gauze swabs and forceps provided.
9. With a clean swab stick, apply a *thin* layer of white soft paraffin (supplied) to puncture area to promote capillary blood aggregation. **Remove excess with a new dry swab stick before skin puncture.**
10. The sampling wand comes pre-fitted with an FBS blade. Lightly holding the shaft and rear finger grip (**A**), push the grip forward to expose the blade (**B**). **Ensure the blade is locked forward before use.**

11. Using the Sampling Wand, perform skin puncture through the amnioscope. Either (i) make a vertical stab, rotate the blade 90° and make a second stab to form a cross or (ii) make a vertical stab, then angle the blade downwards and extend the incision using the front edge of the blade to produce a single incision 4mm long.
12. **Always wait for a good sized (3-4mm) bead to form on the surface before collection**. Never slash at the scalp and avoid multiple small stabs.
13. To collect the blood sample, depress the button on the wand shaft to retract the FBS blade and angle the capillary tube into the bead.
14. Rapidly collect blood sample, filling the tube completely, *without air bubbles*. Collection needs to be expedited to avoid blood clotting.
15. Twist lightly to remove the capillary tube from sampling wand and cap both ends with the covers supplied. Take to the analyser **immediately**. Follow the analyser manufacturer's instructions as regards sample mixing.
16. Always ensure the **clean**, non-fetal exposed end is presented to the analyser. Note: some analysers require a 'clot catcher' (not provided). Please refer to the manufacturer's instructions for guidance on when a clot catcher is needed.
17. Always ensure that the **correct** measuring mode is selected *before* performing analysis.

Post sampling:

18. Haemostasis is accomplished by applying pressure to the puncture site with a swab for one full contraction.

Disposal: This device, its accessories and the consumables used with it, should be handled and disposed of in accordance with policy of the healthcare setting and with regard to all applicable regulations, including but without limitation to, those pertaining to human health and safety, and care of the environment. Ensure sharps are disposed of in sharps container and devices should be decontaminated before disposal or transfer to a third party and supplied with a certificate of decontamination. The Rocket Medical Amnilume™ light module contains lithium-ion based power cells. Lithium cells are suitable for incineration or landfill disposal.

Waste Electrical and Electronic Equipment: In Europe, disposal of this device must be undertaken with regard to the WEEE directive (2002/96/EC). Do not dispose of this device in municipal waste.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

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Copyright © 2023 Alle Rechte vorbehalten.
Copyright © 2023 Todos los derechos reservados.
Copyright © 2023 Kaikki oikeudet pidätetään.
Copyright © 2023 Tous droits réservés.
Copyright © 2023 Alle rechten voorbehouden.