

Lucy and Lucy's Mum- Instrumental Birthing Simulator


Lucy and Lucy's Mum are separately available parts of the Instrumental Birthing Simulator Trainer.

Lucy accurately mimics the average size and shape of the fetal head at birth. The head is complete with articulating cranial plates and a scalp.

Lucy's Mum has a full pelvis, the joints of which are slightly loosened to simulate the movement in those joints during childbirth.

The fetal head must travel around the Curve of Carras and its attitude. The position can be changed to simulate changes in the attitude and position for birth requiring intervention.

The Instrumental Birth Simulator also allows for the practice of forceps deliveries as well as the practicing of vacuum extraction techniques.



Ordering Information

Lucy's Mum	0 T 003
Lucy- Fetal head with Scalp	VAC 6009S

Lucy Pack Contains:

Lucy's Mum training model, Lubricant, Gloves, Talcum powder, Full instructions, Interactive e-training
--

